

Prepare for COVID-19 Testing

Rapid, reliable SARS-CoV-2 testing at the Point-of-Care easy-to-use BD Veritor™ Plus System provides reliable SARS-CoV-2 results in 15 minutes

When your patients are in need of fast, reliable SARS-CoV-2 testing, turn to the BD Veritor™ Plus System. Offering lab-quality results at the point of care, in a simple-to-operate, handheld instrument with results in 15 minutes.

Simplify the testing process

- Offers mid nasal swabbing technique for a less invasive experience
- Easy operation and 1-button functionality may potentially reduce manual test processing errors
- Enables intuitive sample processing with prefilled, unitized tubes color-coded by reagent or assay type

Achieve fast, reliable results


- Displays easy-to-read digital results for SARS-CoV-2 in 15 minutes
- Records results on secured drive
- Advanced particle technology enhances sensitivity by using a proprietary process to produce highly stable, modified colloidal metal particles, helping improve test performance*
- Adaptive read technology helps improve specificity to reduce false-positive results by compensating for background and nonspecific binding

Provide workflow efficiency

- Adapts easily to your workflow by offering 2 operational modes
 - **Walk away:** Test device is inserted immediately into Analyzer, enabling staff to multitask while sample incubates
 - **Analyze now:** Test device is inserted after incubation time is complete, allowing batches of samples to be tested

Clinical Performance¹

PPA	NPA
84%	100%
(95% C.I. 67%–93%)	(95% C.I. 98–100%)



Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management.

Provide result traceability

- Download and display LOT number, Patient/specimen ID, Operator ID and test records with BD Veritor™ InfoScan module
- Offers result-printing capabilities via USB port

*in comparison to visually read tests

For more information, contact your Cardinal Health Sales Representative or email PatientCareSolutions@cardinalhealth.com

The availability of collaborative practice agreements are dependent on state regulations.

© 2020 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and ESSENTIAL TO CARE and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 1R118-812529 (06/2018)


CardinalHealth
Essential to care™

cardinalhealth.com